Nederlands	2
Français.....	6
Deutsch.....	10
English	14

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LAPAROX®

Samenstelling, toedieningsvorm:

Kooldioxide (CO₂), 100%

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Composition, presentation:

Carbon dioxide (CO₂), 100 %

LAPAROX is packaged as liquefied gas under pressure in gas cylinders or as bulk. 1 liter of LAPAROX corresponds to 0,75 kg at a temperature de 15°C.

Intended application:

LAPAROX is intended to be used as an accessory to a CE-marked medical device for administration of carbon dioxide.

Common medical procedures where LAPAROX is used are insufflation for creation of visibility during endoscopic, minimal invasive surgical procedures and also cryosurgery.

Special warnings/possible side-effects:

LAPAROX should always be administered in accordance with the instructions from the medical device manufacturer, by educated personnel.

CO₂ (a gas heavier than air) can accumulate on low points and make the atmosphere dangerous.

Good ventilation is required in areas where LAPAROX is used.

Possible side effects of LAPAROX are carbon dioxide retention, cardiovascular effects, pneumothorax, emphysema, venous gas emboli and hypoxia.

Contra-indications diagnostic Laparoscopy:

- Inability to tolerate pneumoperitoneum or anesthesia
- Patients with an obvious indication for surgical intervention
- Hemodynamic instability (defined by most studies as systolic pressure < 90 mm Hg)
- Patients with an uncorrectable coagulopathy or uncorrectable hypercapnia >50 torr
- Patients with a tense and distended abdomen (i.e., clinically suspected abdominal compartment syndrome)
- Patients with abdominal wall infection
- Known or obvious intra-abdominal injury
- Patients with extensive previous abdominal surgery with multiple incisional scars or after a laparotomy within the last 30 days
- Multiple adhesions
- For diagnostic staging laparoscopy:
 - In pancreatic adenocarcinoma or esophageal tumors: Known metastatic disease
 - In colorectal cancer: Patients with resectable liver metastases from colorectal cancer but with no evidence of extrahepatic disease on non-invasive imaging
 - In gastric cancers: Patients with early stage gastric cancer (T1 or T2) or gastric cancers complicated by obstruction (palliative surgery is needed)
 - For primary hepatic tumors: Patients with known unresectable hepatic disease such as major vessel or organ invasion, which are not candidates for surgery
 - For biliary tract tumors: known metastatic or unresectable disease, known stage T1 disease found incidentally may potentially be treated with cholecystectomy alone

Contra-indications cryosurgery:**Absolute:**

- Blood dyscrasias of unknown origin.
- Cold intolerance.
- Raynaud's disease.
- Cold urticaria.
- Cryoglobulinemia.
- Lesions in which tissue pathology is required.
- Lesions in areas of compromised circulation.
- Sclerosing basal cell carcinoma (BCC) or recurrent BCC or squamous cell carcinoma (SCC) located in high-risk areas like the temple or nasolabial folds.

Relative:

- Keloidal tendency.
- Collagen vascular diseases.
- Dark-skinned individuals due to the high risk of developing cosmetically unacceptable protracted hypopigmentation.
- Lesions over the nasolabial fold, eyelid margins, ala nasi and hair-bearing areas (high risk of developing alopecia, especially cicatricial alopecia).
- Patients with sensory loss at lesional sites.
- Pyoderma gangrenosum.

Precautions for use:

Contains (liquified) gas under pressure; may explode if heated.

Handle material with clean hands.

Do not transfer pressurized gas from one cylinder to another.

Do not stand in front of the valve outlet or behind the pressurized gas cylinder, but always at the side.

Use standard and specific connection fittings or hoses for medical CO₂.

Unless explicitly specified, use the cylinder in an upright position.

The use of the product at high flow can cause icing at the cylinder and connections.

Do not attempt to repair a defective valve.

Do not disconnect or try to tighten the connection without first depressurizing the system.

Instruction for use:

Before use, refer to the manufacturer's instructions of the device associated with the medical CO₂ and follow the recommendations of the healthcare professional.

Cylinders:

At delivery and prior to any use, check the package integrity and the tamper evidence sealing system.

In particular, check whether the equipment is suitable and in good condition. Ensure that the residual amount will allow continuity of supply throughout the product delivery process by weighing the cylinder. Since LAPAROX is a liquid gas under pressure, the pressure drops quickly as soon as the cylinder no longer contains any liquid.

If necessary, clean and decontaminate the cylinder and check its integrity.

Remove the seal (tamper evidence) from the package, check that the valve is clean and that no foreign particles are present.

Open the cylinder valve slightly until the gas flows out and immediately close the valve again (this will remove any particles). Repeat this step 2 to 3 times.

Connect the cylinder to the corresponding device:

- Screw connection: align the regulator and cylinder connections.
- Pin index connection: align the locating pins of the valve to the holes.

Secure the connection by hand until it reaches the stop.

Before opening the valve check that the couplings are properly fitted.

Do not open the cylinder valve if not connected (risk of cold burns).

Slowly open the cylinder valve completely, without forcing.

Suddenly opening the valve may lead to liquifying the product and could pose a risk of cryogenic burning.

After use, always close the cylinder valve and depressurizing the system.

Bulk:

Slowly open the valve completely, without forcing.

Suddenly opening the valve may lead to liquifying the product and could pose a risk of cryogenic burning.

After use, always close the cylinder valve and, if applicable, depressurise the system.

Special recommendations:

For medical use by qualified personnel. Use should be by competent and qualified personnel according to the instructions for use with the administration device.

When transporting the cylinders of LAPAROX, secure the cylinders to avoid any risk of shocks or falling.

Handle cylinders with caution and particularly those with a capacity of <5 liters (no protective cap on Pin-Index taps).

Do not lift the cylinder by its valve. Keep empty cylinders with valve closed.

Storage/transport:

Handle carefully.

Do not exceed the expiry date which is indicated on the package label.

Keep in well-ventilated areas, away from sources of heat, bad weather and at temperatures below 50°C.

Storage areas must be sufficiently ventilated.

In storage, during transport and in use: secure the cylinders with an appropriate means (chains, hooks, ...) to keep them in a vertical position and avoid any risk of accidental falling over.

Store and transport upright with closed valve and pressure regulator removed.
Not MRI compatible.

Possible incidents:

Gas leak: close the valve, purge the coupling system. If the leak persists, purge the system from the cylinder, allow the gas to be evacuated without attempting to seal or repair the valve and return the defective packaging to the manufacturer.

Adverse event must be reported to the national competent authority in accordance with the national vigilance system for medical devices:

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Further Information:

This instructions for use was last updated: January 2021.

Please refer for further information relating to cryogenic gases and Liquid Carbon Dioxide to the MSDS available from www.linde-healthcare.nl/shop or contact the manufacturer (see below).

Manufacturer:

Linde Gas Therapeutics Benelux B.V.
De Keten 7
5651 GJ Eindhoven
The Netherlands

The Netherlands:	Tel. +31 (0)40 2825825
Belgium/Luxembourg:	Tel. +32 (0)70 233826